

Nemera

rommelag®

Novelia® BFS

THE BEST OF TWO WORLDS:

PRESERVATIVE-FREE MULTIDOSE NOVELIA® COMBINED WITH BOTTELPACK®-TECHNOLOGY



Multidose ophthalmic closing tip system avoids the need of preservatives and prevents bacterial contamination during use. Advanced Blow-Fill-Seal (BFS) ensures efficient and aseptic processing.

Key benefits of Novelia®

- **Preservative-free** to respect the ocular surface
- **100% controlled and safe** thanks to its patented **PureFlow™ Technology**
- Works with high viscosity suspensions and solutions up to 1500cps
- **User-friendly** and **intuitive**: easy eye targeting **blue tip** feature, as easy to use as any standard eyedropper
- **One drop at a time** into the patient's eye
- **Calibrated drops**, precise dose for a better adherence to the medication



Key benefits of bottelpack® BFS


- **Compact** BFS processing lines with **highest aseptic standards**
- Proven and **cost-effective** production
- Thousands of bottelpack® machines sold worldwide for over 50 years
- **Customized containers** for differentiated products
- **Easy to squeeze** bottles



Nemera

rommelag®

FUNCTIONING OF THE NOVELIA® SYSTEM



THE BOTTELPACK® BLOW-FILL-INSERT-SEAL PROCESS



Novelia® BFS prototypes are available.
Please contact us for more information.



www.nemera.net

20, Avenue de la Gare - 38290 LA VERPILLIERE - FRANCE
Tel: +33 (0)4 74 94 06 54 - Fax: +33 (0)4 74 94 90 60

Contact : fanny.sellier@nemera.net

www.rommelag.com

Fabrikweg 16 - CH-5033 Buchs - SWITZERLAND
Tel: +41 62 834 55 55 - Fax: +41 62 834 55 00

Contact : boris.kovacic@rommelag.ch

we put
patients
First