

Intimate care & Oral administration



- VAGINAL AND RECTAL APPLICATORS
- ORAL DOSERS, DROPPERS AND DRIPPERS



Nemera established an operational footprint in Brazil in 2021 and expanded its products and services to better serve the pharmaceutical industry in Latin America.

Our devices aim to be used for 2 main therapeutic areas:

INTIMATE CARE:

- Vaginal applicators
- Rectal applicators

ORAL ADMINISTRATION:

- Oral dosers
- Oral droppers
- Oral drippers





Intimate Care

Gynecological pathologies with highest prevalence [13-50 years old]

VAGINAL CANDIDIASIS¹

- Prevalence: ≈60% women
- Fungal infection causing irritation, discharge and itchiness of vagina and vulva

BACTERIAL VAGINOSIS¹

- Prevalence: ≈30% women
- Vaginal infection caused by excessive growth of bacteria. Common symptoms include increased vaginal discharge and unpleasant vaginal odor



Vaginal atrophy during menopause [45-60 years old]

VAGINAL ATROPHY²

- Estrogen deficiency occurring after menopause accelerates the deterioration of vulvovaginal tissue
- Vaginal dryness, inflammation, dyspareunia

Hemorrhoids and anal fissures

- Prevalence: ≈ 4.4% in the general population
- Increases with age, with a peak in persons aged 45-65 years

1) High Global Burden and Costs of Bacterial Vaginosis: A Systematic Review and Meta-Analysis. Peebles, Kathryn MPH*; Vellozo, Jennifer MPH*; Balkus, Jennifer E. PhD*†§; McClelland, R. Scott MD*†§; Barnabas, Ruanne V. Sexually Transmitted Diseases 46(5):p 304-311, May 2019 |

2) Fabion SS, Sood R, Kapoor E. Genitourinary Syndrome of Menopause: Management Strategies for the Clinician. Mayo Clin Proc. 2017 Dec;92(12):1842-1849

Our products portfolio

VAGINAL APPLICATORS FOR LIQUID (Gel, Cream)

- 26 different locks making it possible to offer dosages between 0,5 to 6g
- Available in different colors and models
- Individually packaged in cartons or sachets, according to the need for each drug



IN DEVELOPMENT

VAGINAL APPLICATORS FOR SOLID (Pills, Tablets)

RECTAL APPLICATORS

- Packaged in cartons, or sachets, according to the need for each drug
- Available in different colors



Oral administration

Pediatric pathologies & treatments

[baby & child]

ANALGESIC NON-NARCOTIC

- Fever
- Pain

ANTI-ACID & ANTI-FLATULENT

- Reflux
- Infant colic pain



Medical food supplementation



LIQUID FOOD SUPPLEMENT

- Vitamins (D, C, B...)
- Minerals (calcium, iron, magnesium...)
- Botanicals (ginger, medical cannabidiol...)
- Omega fatty acids

Our products portfolio

ORAL DOSERS

- Available for 3, 5, 10 and 15ml
- Individually packaged
- Availability of covers and adapters in sizes, 11, 13, 16, 19, 22 and 26mm



ORAL DROPPERS

- Available in different colors
- With or without protective cover
- Can be packaged individually
- 51, 60 and 72mm cannula, graduated and not graduated

ORAL DRIPPERS

- Available in different colors and models
- Size for 18mm termination bottles



Others

CAPS 18MM

- Available in different colors
- Different seal options
- Also available in version with hole, for droppers



SPRAY CAPS

- Produced with resistant material to surface treatments

Our services

ANALYTICAL SERVICES & DESIGN VERIFICATION

- Fit for purpose testing plan and set-up development
- Testing results formatted
- Device verification testing report



REGULATORY & COMPLIANCE DOSSIER SUPPORT

- Strategy, planning, classification, regulations
- Standards, guidance and pharmacopeia interpretation and analysis
- Risk management and reliability
- Clinical evaluation
- Regulatory registration and CE marking
- Notified body opinions
- Post market activities (plan, report...)



Learn more about
our services at:
nemera.net/our-services

Our broad offering

Offering a broad range of solutions in major delivery routes, we are the holistic development and manufacturing partner for your custom drug delivery device solution.

